



SPECIFICATION SHEET

Product Name : Acetonitrile with 0.1% Trifluoroacetic acid, LC-MS
Product Number : 901.088
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Assay (GC)	≥ 99.5 %	%
Assay (Trifluoroacetic acid)	0.093 - 0.107 %	%
Appearance	Clear colorless liquid	
Calcium (Ca)	≤ 0.00005 %	%
Potassium (K)	≤ 0.00005 %	%
Magnesium (Mg)	≤ 0.00005 %	%
Sodium (Na)	≤ 0.0002 %	%
Transmittance at 230 nm	≥ 50 %	%
Transmittance at 254 nm	≥ 90 %	%
LC-MS application	passes test	
Grad. elution H.Peak at 254 nm	≤ 0.05 AU	AU
Filtered through 0.2 μm		

Quality Management
Isolab Laborgerate GmbH

