



SPECIFICATION SHEET

Product Name : Ammonium acetate Extra pure
Product Number : 903.022
Brand : ISOLAB
CAS Number : 631-61-8
Formula : $\text{CH}_3\text{COONH}_4$
Formula Weight : 77.08 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay	≥ 96.0	%
Identity	passes test	
pH-value (5 %; water, 25 °C)	6.7 - 7.3	
Chloride (Cl)	≤ 0.002	%
Sulphate (SO_4)	≤ 0.01	%
Heavy metals (as Pb)	≤ 0.0005	%
Fe (Iron)	≤ 0.001	%
Residue on ignition (as sulphate)	≤ 0.02	%
Water	≤ 2.5	%

Quality Management
Isolab Laborgerate GmbH

