



SPECIFICATION SHEET

Product Name : Butylated Hydroxyanisole
Product Number : 908.069
Brand : ISOLAB
CAS Number : 25013-16-5
Formula : $C_{11}H_{16}O_2$
Formula Weight : 180.25 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	White to light yellow, powder	
Solubility	Freely soluble in alcohol, propylene glycol. Insoluble in water	
Assay (HPLC)	≥ 98.5 %	%

Quality Management
Isolab Laborgerate GmbH

