



CERTIFICATE OF ANALYSIS

Product Name : Buffer Solution pH 10.00 ± 0.01 at 20°C
(Di-sodium Tetraborate, 1M Sodium Hydroxide)
Product Number : 908.B10
Batch Number : LR2453511AGW
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -
Production Date : 11/2017
Expiration Date : 11/2019

TEST NAME	SPECIFICATION	UNIT	RESULTS
pH	9.99 - 10.01 at 20°C	-	10.00 at 20°C
Measured with a combination glass electrode after multiple point calibration with reference materials. It is certified traceable to the following NIST(USA), SRM 191d-I Sodium Bicarbonate, SRM 191d-II Sodium Carbonate, SRM 186-I-g Potassium Dihydrogen Phosphate and SRM 186-II-g Disodium Hydrogen Phosphate.			

Quality Management
Isolab Laborgerate GmbH

