



SPECIFICATION SHEET

Product Name : Ethyl acetate, LC-MS
Product Number : 920.018
Brand : ISOLAB
CAS Number : 141-78-6
Formula : C₄H₈O₂
Formula Weight : 88.11 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (GC)	≥ 99.7 %	%
Non-volatile matter	≤ 0.0005 %	%
Water (Karl Fischer)	≤ 0.05 %	%
Free alkali (as NH ₃)	≤ 0.0005 %	%
Calcium (Ca)	≤ 0.00001 %	%
Potassium (K)	≤ 0.00001 %	%
Magnesium (Mg)	≤ 0.00001 %	%
Sodium (Na)	≤ 0.00001 %	%
Transmittance at 260 nm	≥ 50 %	%
Transmittance at 300 nm	≥ 98 %	%
Suitable for HPLC-MS	passes test	

Quality Management
Isolab Laborgerate GmbH

