



### SPECIFICATION SHEET

**Product Name** : Ethylene glycol, Extra pure  
**Product Number** : 920.032  
**Brand** : ISOLAB  
**CAS Number** : 107-21-1  
**Formula** :  $C_2H_6O_2$   
**Formula Weight** : 62.07 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity	≥ 99 %	%
Appearance	Colorless, liquid	
Infrared spectrum	passes test	
Color test	≤ 10 APHA	
Acidity	≤ 0.1 ml N %	%
Sulphated Ash	≤ 0.01 %	%

Quality Management  
Isolab Laborgerate GmbH

