



SPECIFICATION SHEET

Product Name : Formic acid 85%, Extra pure
Product Number : 923.022
Brand : ISOLAB
CAS Number : 64-18-6
Formula : CH_2O_2
Formula Weight : 46.03 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Colorless, clear, liquid	
Assay	$\geq 85.0 \%$	%
Color (PT-CO)	≤ 10	
Cl (Chloride)	$\leq 0.002 \%$	%
SO ₄ (Sulfate)	$\leq 0.001 \%$	%
Fe (Iron)	$\leq 0.0001 \%$	%
Evaporation residue	$\leq 0.006 \%$	%

Quality Management
Isolab Laborgerate GmbH

