



SPECIFICATION SHEET

Product Name : n-Heptane, Extra pure
Product Number : 930.012
Brand : ISOLAB
CAS Number : 142-82-5
Formula : C₇H₁₆
Formula Weight : 100.2 g/mol

| TEST NAME | SPECIFICATION | UNIT |
|--|---------------|-------|
| Purity (GC) | ≥ 99.0 % | % |
| Identity (IR) | passes test | |
| Free acid (as CH ₃ COOH) | ≤ 0.005 % | % |
| Density (d 20°C/ 4°C) | 0.683 - 0.684 | |
| Matter discoloured by H ₂ SO ₄ | ≤ 10 Hazen | Hazen |
| Evaporation residue | ≤ 0.005 % | % |

Quality Management
Isolab Laborgerate GmbH

