



SPECIFICATION SHEET

Product Name : Hydrogen peroxide 35 %, Extra pure
Product Number : 932.122
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Concentration (± 0.2)	35 %	%
Appearance	Clear colourless liquid	
Density (g/cm^3) 25 °C	1.125 - 1.130	g/cm^3
Color (APHA)	< 10	
Conductivity	> 50 $\mu\text{S/sec}$	$\mu\text{S/sec}$
Acidity	≤ 0.03 %	%
Stability	≤ 0.3 ml/min	ml/min
Chloride (Cl ⁻)	≤ 0.0001 %	%
Nitrate (NO ₃)	≤ 0.015 %	%
Phosphate (PO ₄)	≤ 0.005 %	%
Sulfate (SO ₄)	≤ 0.001 %	%
As (Arsenic)	≤ 0.0003 %	%
Cu (Copper)	≤ 0.0001 %	%
Fe (Iron)	≤ 0.00005 %	%
Mn (Manganese)	≤ 0.0001 %	%
Pb (Lead)	≤ 0.001 %	%
Sn (Tin)	≤ 0.001 %	%
Non volatile matter	≤ 0.015 %	%

Quality Management
Isolab Laborgerate GmbH

