



### SPECIFICATION SHEET

**Product Name** : May and Grunwald's stain, for microscopy  
**Product Number** : 946.D02  
**Brand** : ISOLAB  
**CAS Number** : -  
**Formula** : -  
**Formula Weight** : -

TEST NAME	SPECIFICATION	UNIT
Appearance	A dark green powder	
Loss on drying (at 105°C)	≤ 10 %	%
Absorbance maxima at λ1	648-653 nm	nm
Absorbance maxima at λ2	520 - 525 nm	nm
Specific absorptivity A 1%/1 cm at λ1	1100 - 1300	
Specific absorptivity A 1%/1 cm at λ2	600 - 800	

Quality Management  
Isolab Laborgerate GmbH

