



SPECIFICATION SHEET

Product Name : Oil red O, (C.I. 26125) for electrophoresis
Product Number : 955.D01
Brand : ISOLAB
CAS Number : 1320-06-5
Formula : $C_{26}H_{24}N_4O$
Formula Weight : 408.5 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Brownish red powder	
Absorption maxima(λ) (Chloroform)	528 - 534 nm	nm
Loss on drying 105°C	≤ 5 %	%

Quality Management
Isolab Laborgerate GmbH

