



SPECIFICATION SHEET

Product Name : n-Pentane, Extra pure
Product Number : 959.012
Brand : ISOLAB
CAS Number : 109-66-0
Formula : C₅H₁₂
Formula Weight : 72.15 g/mol

| TEST NAME | SPECIFICATION | UNIT |
|-------------------------------------|---------------|------|
| Purity (GC) | ≥ 99.0 % | % |
| Identity (IR) | passes test | |
| Free acid (as CH ₃ COOH) | ≤ 0.001 % | % |
| Non volatile matter | ≤ 0.005 % | % |
| Readily carbonizable substances | passes test | |

Quality Management
Isolab Laborgerate GmbH

