

SPECIFICATION SHEET

Product Name : Phloxin B, (C.I.45410) for microscopy
Product Number : 959.D01
Brand : ISOLAB
CAS Number : 18472-84-2
Formula : $C_{20}H_2Br_4Cl_4Na_2O_5$
Formula Weight : 829.64 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Red coloured powder	
Absorption maxima in 50% ethanol	546 - 550 nm	nm
Absorptivity (1cm, λ_{max})	≥ 750	
Loss on drying (at 110°C)	$\leq 10 \%$	%

Quality Management
 Isolab Laborgerate GmbH


