



### SPECIFICATION SHEET

**Product Name** : Potassium hexacyanoferrate(II) trihydrate, for analysis ACS,ISO,Reag. Ph Eur  
**Product Number** : 960.076  
**Brand** : ISOLAB  
**CAS Number** : 14459-95-1  
**Formula** :  $C_6FeK_4N_6 \cdot 3H_2O$   
**Formula Weight** : 422.39 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (permanganometric)	99 - 102 %	%
Insoluble matter	$\leq 0.005$ %	%
Chloride (Cl)	$\leq 0.01$ %	%
Sulphate (SO <sub>4</sub> )	$\leq 0.005$ %	%
Cd (Cadmium)	$\leq 0.0005$ %	%
Cu (Copper)	$\leq 0.002$ %	%
Pb (Lead)	$\leq 0.002$ %	%
Na (Sodium)	$\leq 0.01$ %	%

Quality Management  
Isolab Laborgerate GmbH

