



SPECIFICATION SHEET

Product Name : Polysorbate 20, For synthesis
Product Number : 960.151
Brand : ISOLAB
CAS Number : 9005-64-5
Formula : $C_{58}H_{114}O_{26}$
Formula Weight : 1228 g/mol

TEST NAME	SPECIFICATION	UNIT
Density (d 20°C/ 4°C)	1.095 - 1.105	
Saponification value	40 - 50	
Hydroxyl value	96 - 108	
Identity (IR)	passes test	

Due to its specific melting range the product may be solid, liquid, a solidified melt or a supercooled melt.

Quality Management
Isolab Laborgerate GmbH

