



SPECIFICATION SHEET

Product Name : Sodium Disulfite Ph Eur,BP,JP,NF,E 223
Product Number : 969.073
Brand : ISOLAB
CAS Number : 7681-57-4
Formula : $\text{Na}_2\text{O}_5\text{S}_2$
Formula Weight : 190.11 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay ($\text{Na}_2\text{S}_2\text{O}_5$)	≥ 98.0 - 99.5	%
Assay (SO_2)	≥ 66.0 – 67.0	%
Identity	passes test	
Appearance of solution	passes test	
pH (10% sol)	4.0 - 4.5	
Chloride (Cl)	≤ 0.005	%
Heavy metals (as Pb)	≤ 0.001	%
As (Arsenic)	≤ 0.00005	%
Fe (Iron)	≤ 0.0005	%
Hg (Mercury)	≤ 0.00001	%
Pb (Lead)	≤ 0.00005	%
Se (Selenium)	≤ 0.00005	%

Quality Management
Isolab Laborgerate GmbH

