



SPECIFICATION SHEET

Product Name	:	Sodium hydroxide, pellets pure
Product Number	:	969.112
Brand	:	ISOLAB
CAS Number	:	1310-73-2
Formula	:	NaOH
Formula Weight	:	40.00 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (acidimetric, NaOH)	≥ 99 %	%
Carbonate (as Na ₂ CO ₃)	≤ 1.0 %	%
Chloride (Cl)	≤ 0.005 %	%
Nitrogen Compounds(N)	≤ 0.0005 %	%
Phosphate (PO ₄)	≤ 0.02 %	%
Silicate (SiO ₂)	≤ 0.01 %	%
Sulphate (SO ₄)	≤ 0.005 %	%
Al (Aluminium)	≤ 0.002 %	%
Ca (Calcium)	≤ 0.005 %	%
Copper (Cu)	≤ 0.0005 %	%
Fe (Iron)	≤ 0.001 %	%
Ni (Nickel)	≤ 0.001 %	%
K (Potassium)	≤ 0.05 %	%
Heavy Metals (as Pb)	≤ 0.0005 %	%

Quality Management
Isolab Laborgerate GmbH

