



SPECIFICATION SHEET

Product Name : Sodium hydroxide 0.01 mol/l (0.01 N) , Volumetric solution  
Product Number : 969.19V  
Brand : ISOLAB  
CAS Number : -  
Formula : -  
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Amount of substance concentration	0.00995 - 0.01005 mol/l	mol/l
Traceability	NIST	
Indicator	phenolphthalein	

Quality Management  
Isolab Laborgerate GmbH

